

Designed for the Iliofemoral Veins

The Venovo™ Venous Stent and Atlas™ Gold PTA Dilatation Catheter are purpose-bulit for the veins with the sizes and strengths needed to treat iliofemoral venous lesions.

Size Ranges

Diameter (mm)	Length (mm)							
	20	40	60	80	100	120	140	160
10		●	●	●	●	●	●	●
12	●	● ●	● ●	●	●	●	●	●
14	●	● ●	● ●	●	●	●	●	●
16	●	● ●	● ●	●	●	●	●	●
18	●	● ●	● ●	●	●	●	●	●
20	●	● ●	●	●	●	●	●	●
22	●	●						
24	●	●						
26	●	●						

Atlas™ Gold
PTA Dilatation Catheter

Venovo™
Venous Stent System

NOTE: The ATLAS™ GOLD PTA Dilatation Catheter has rated burst pressures up to 18 mm. The VENOVO™ Venous Stent System is indicated for the treatment of symptomatic iliofemoral venous outflow obstruction. The ATLAS™ GOLD PTA Dilatation Catheter is indicated for use in percutaneous transluminal angioplasty of the peripheral vasculature, including the iliac arteries and iliac and femoral veins, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for post-dilatation of stents and stent grafts in the peripheral vasculature. This catheter is not for use in coronary arteries.

Please consult product labels and instructions for use for indications, contraindications, hazards, warnings, and precautions.

Please consult product labels and instructions for use for all indications, contraindications, hazards, warnings and precautions.

바드코리아(주) 06236 서울시 강남구 테헤란로 142, 16층
T. 02-2188-2900 F. 02-539-3081

crbard.com

